

Chronomask®

Online sale on innoralis.fr

And the Class III treatments become a child's play!



› Light Chronomask®

› A COMFORTABLE AND STREAMLINED DESIGN

Anatomic chin rest, forehead rest completely flexible, field of vision completely free and metallic frame reduced at its minimum make this orthopaedic mask the lightest and most comfortable on the market.

› UNIQUE SIZE EASY TO ADJUST

Its measurements have been studied in order to be convenient as much for the young children's as the teenagers. A forehead Velcro, an adaptable frame and a discrete set screw of the traction height make the Chronomask® an orthopaedic mask quickly adjustable with an easy managing stock.

› WASHABLE HYGIENIC REST PAD

Forehead and chin pads are made with a soft and comfortable textile, which is washable in washing machine and interchangeable. The Chronomask® is delivered with a restock of textile pads to maintain in the easiest way these parts which are in contact with the skin of our young patients.



keep reading ›

➤ Connected Chronomask®



› RECORDING OF THE WEARING DATAS

A sensor records the wearing time periods of the mask and the intensity of the used forces.

› PLAYFUL AND MOTIVATING SMARTPHONE APP FOR PATIENTS

Each week the connected Chronomask® Bluetooth module send its datas to a motivating smartphone / tablet application belonging to the young patient. If the wearing time is over 12h, the children will move forward through 36 adventures on the Ekyoki tracks... a young panda whom need the energy coming from the Chronomask® wearing time, to find his parents again.

› MONITORING SOFTWARE FOR PRACTITIONERS

At each appointment, the clinician can control the correct use of his medical prescriptions by using a software that presents, under average forms and graphical representations, the wearing datas of the mask since the previous appointment.



“ These technological advances offer to the clinician a total control over the orthopaedic treatment : say goodbye to wasted time caused by relapses, and welcome to modern orthopaedic and efficiency ! ”

Dr. Jean-Louis Raymond,
Chronomask® Inventor.



innOralis®

Le Millenium, 18 rue Albert Camus
90000 Belfort – France
contact@innoralis.fr



Join the innovation community at **innOralis®.fr**